

EvaDrive Dimensions

| Name | EvaDrive | |
|--|--|-------------------------------------|
| Part Number | 60900002 | 60900003 |
| Description | EvaDrive450 | EvaDrive600 |
| Unit | 1 piece | 1 piece |
| Base width, outer dimensions | 27in – 38.6in (686mm – 980mm) | 29.5in – 45.2in (749mm – 1148mm) |
| Base width, inner dimensions | 23in – 34.6in (584mm – 879mm) | 25.6in – 41.3in (650mm – 1049mm) |
| Sling bar width | 17.7in (450mm) | 23.6in (600mm) |
| Lifting speed | 1.4in/sec (36mm/sec) | |
| Driving speed | 4.3mph (7km) | |
| Castors | 3.9in (99.0mm) (front) 8in (203.2mm) (rear) | |
| Safe Working Load | 450lb (204kg) | 600lb (272kg) |
| Overall weight of lift (including batteries) | 128lb (58kg) | 141.1lb (64kg) |