

# WE'RE THE COMPANY THAT BRINGS YOU KALEI...



## AND KALEI BRINGS YOU THE P-440

Unique and innovative, the stylish Prism Medical P-440 stretches the boundaries of portable ceiling lifts. Ceiling lifts dramatically reduce injuries resulting from the handling of patients/residents in both institutional and homecare settings. The P-440 has been designed with a full set of standard features and can be operated safely by a single caregiver.


The P-440 is designed to be light, easy to hold and is capable of lifting up to 440 lbs (200 kgs).



10888 Metro Court, St. Louis, MO 63043 | 866.891.6502 | [www.handicareusa.com](http://www.handicareusa.com)  
485 Millway Avenue, Unit 2, Concord, ON L4K 3V4 | 877.304.5438 | [www.handicare.ca](http://www.handicare.ca)

# P-440

## PORTABLE



## WE'RE THE COMPANY THAT BRINGS YOU PORTABLE FLEXIBILITY.

The P-440 solves many of the problems that make patient transfer complex or difficult. It's lightweight, portable, and designed to provide complete control for the caregiver. In addition, it offers a wide range of safety mechanisms to assure confidence throughout the transfer.

### PRODUCT HIGHLIGHTS

- Lightweight and portable
- Integral full-width carry bar with unique contoured shape for optimal client comfort
- On-board controls allow the caregiver to operate the system. An optional hand control is also available
- The P-440 is charged with a standard plug and socket
- Safety mechanisms include an over-speed governor, built-in safety limits, slack strap sensor and emergency lowering
- Patient safety is assured with positive locking latches, securing the sling each and every time
- The P-440 offers a comfortable built-in ergonomic handle, allowing the caregiver to carry and transport with ease
- The handle doubles as a connection point to allow for room to room transfers eliminating the need for costly home modifications

### TECHNICAL SPECIFICATIONS

- Charger Output: 24 VDC, 1.5 Amps
- Batteries: 24 VDC (2 x 12 VDC), 2.3 AH
- Lift Case: Flame Retardant ABS
- Hand Control: Pneumatic
- Lifting Range: Up to 80" (200 cm)
- Unit Weight: 10.5 lbs (4.8 kgs)
- Maximum Load: 440 lbs (200 kgs)



10888 Metro Court, St. Louis, MO 63043 | 866.891.6502 | [www.handicareusa.com](http://www.handicareusa.com)  
485 Millway Avenue, Unit 2, Concord, ON L4K 3V4 | 877.304.5438 | [www.handicare.ca](http://www.handicare.ca)



For Your Convenience

