

ReTurnBelt Dimensions

Name	ReTurnBelt*, nylon with 4 handles							
Part Number	6031	6032	6033	6034	6035	6036	6037	6038
Description	XXS	XS	S	M	L	XL	XXL	XXXL
Unit	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
Dimensions	25.6" – 31.5"	31.5" – 35.4"	35.4" – 41.3"	41.3" – 49.2"	49.2" – 57.1"	57.1" – 61"	61" – 68.9"	68.9" – 76.8"
Safe Working Load	450lbs	450lbs	450lbs	450lbs	450lbs	450lbs	450lbs	450lbs

*Patent: SE-0401517-8

Name	ReTurnBelt*, disposable with 2 handles				ReTurnBelt*, wipeable with 2 handles			
Part Number	6233	6234	6235	6236	6133	6134	6135	6136
Description	S	M	L	XL	S	M	L	XL
Unit	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
Dimensions	35.4" – 41.3"	41.3" – 49.2"	49.2" – 57.1"	57.1" – 61"	35.4" – 41.3"	41.3" – 49.2"	49.2" – 57.1"	57.1" – 61"
Safe Working Load	450lbs	450lbs	450lbs	450lbs	450lbs	450lbs	450lbs	450lbs

*Patent: SE-0401517-8