



SINA Comfort Shower Trolley

A Comfortable Simplified Shower Experience



The Solution for **Safe and Comfortable Showering**

The SINA Comfort Shower Trolley provides residents with a safe, comfortable and dignified bathing experience. The transfer from bed to shower trolley is simplified by moving the SINA next to the bed and using a transfer aid to move the resident onto the shower trolley. The SINA Comfort comes standard with electric height adjustment, enabling easy adjustment to the working height of the nurse for simplified care. The SINA can also be used as a mobile changing table due to its scissor-system which provides stable support to the bed surface right up to the edge. Developed with care homes in mind, the SINA Comfort can fit into any standard bathroom and its overall appearance provides a modern aesthetic.

SINA Comfort is available with a safe working load of 165 kg (364 lb) and 250 kg (551 lb) for bariatric residents with low mobility. Also available with hydraulic height adjustment.

Caregiver and Resident Benefits



Flexible water drainage hose



Padded, comfortable shower mattress



Adjustable rail with safety joints for extension of the shower surface



Electric or hydraulic height adjustment (optional)



Castors enables easy maneuverability



Tiltable siderails with maintenance-free safety joints for a trouble-free transfer.



Positive and negative tilt adjustment of the bed surface.



SINA Comfort XXL. Extra-wide and extremely stable with electrical height adjustment for bariatric patients up to 250 kg (550 lb).

SINA Comfort Shower Trolley

Specifications

Part Numbers

Safe Working Load (SWL)

SINA Comfort EL 160 / EL 200

165 kg (364 lb)

SINA Comfort XXL

250 kg (551 lb)

Dimensions (L x W x H)

SINA Comfort EL 160

Length: 1740 mm (68.5 in)

Width: 760 mm (29.9 in)

Height: 540 – 865 mm (21.3 – 34.0 in)

SINA Comfort EL 200

Length: 2040 mm (80.3 in)

Width: 760 mm (29.9 in)

Height: 540 – 865 mm (21.3 – 34.0 in)

SINA Comfort XXL

Length: 1900 – 2040 mm (74.8 – 80.3 in)

Width: 800 – 900 mm (31.5 – 35.4 in)

Height: 720 – 1040 mm (28.4 – 40.9 in)

Castors

100 mm (3.9 in)

Power Supply

Removable 24V battery

Protection Class

IPX4

Product Weight Without Packaging

110 kg (243 lb)

For additional technical data, please refer to the user manual.

930507408: SINA Comfort EL 160

Electrical height adjustment, with dual control box, including mattress

930507308: SINA Comfort EL 200

Electrical height adjustment, with dual control box, including mattress

930510008: SINA Comfort XXL

Electrical height adjustment, with dual control box, including mattress

Accessories / Options

Shower Mattress 200

Description	Part No.
Length: 1950 – 2050 mm (76.8 – 80.7 in), Width: 750 mm (29.5 in), Height: 200 mm (7.9 in)	930508300

Shower Mattress 160

Description	Part No.
Length: 1600 – 1700 mm (74.8 – 80.3 in), Width: 750 mm (29.5 in), Height: 200 mm (7.9 in)	930508310

Head Part Shower Trolley

Description	Part No.
The height adjustment of the head part occurs through a gas spring for continuous adjustment and setting of the head part for more comfort of the patient or resident.	930507500

Multicare-Upholstery (with head part)

Description	Part No.
Height: 40 mm (1.8 in), with quilted head part. Head and foot end with loops. Colour: blue	930508400

Head Wedge for Shower Trolley and Multicare

Description	Part No.
Length: 410 mm (16.1 in), Width: 540 mm (21.3 in), Height: 150 mm (5.9 in), water-resistant welded PVC cover with inside foam-core. Colour: light blue	930508320

HANDICARE

485 Millway Avenue, Unit 2, Concord, ON, L4K 3V4

Toll Free: 877.304.5438 | Fax: 416.260.0256

customerservice.canada@handicare.com

www.handicare.ca

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

customerservice.northamerica@handicare.com

www.handicareusa.com

Are you ready to invest in innovation
and simplify bathing for your residents?

To learn more, call 1-888-637-8155 or
visit www.handicare.ca / handicareusa.com