



Clip Slings

Aiding in Safe, Comfortable, Dignified Transfers

Handicare by Beka offers a wide range of Clip Slings with head support for four-point clip style spreader bars. Our Clip Slings provide support for everyday patient handling, padding for the leg area allows for increased comfort and prevents material roping and creasing under the legs. Transfers to and from beds, wheelchairs, floors, and bathrooms are made effortless and more comfortable for the resident and caregiver.

Bathing Clip Slings are made from a mesh nylon material that easily allow water to pass through while providing a very high level of comfort in key pressure areas such as the legs and shoulder, the perfect choice for showering and bathing. Padded leg areas and head support add to the overall comfort and support for the resident.

Five different sizes are available from Small to XXL.

Clip Slings Part Numbers



Transfer Slings with Loops

Made of nylon material. With head support and padded leg area. Washable at max 60° C (140° F).

Part No.	Size	Weight	Colour
922005103	S	35–65 kg (77–143 lb)	Red
922005203	M	55–85 kg (121–187 lb)	Yellow
922005303	L	75–125 kg (165–275 lb)	Green
922005403	XL	120–160 kg (264–352 lb)	Blue
922005503	XXL	150–230 kg (330–506 lb)	Blue



Bathing Slings with Loops

Made of mesh nylon material. With head support and padded leg area. Washable at max 60° C (140° F).

Part No.	Size	Weight	Colour
922006103	S	35–65 kg (77–143 lb)	Red
922006203	M	55–85 kg (121–187 lb)	Yellow
922006303	L	75–125 kg (165–275 lb)	Green
922006403	XL	120–160 kg (264–352 lb)	Blue
922006503	XXL	150–230 kg (330–506 lb)	Blue



Toilet Slings with Loops

Made of nylon material. With head support and padded armpit and leg area. Washable at max 60° C (140° F).

Part No.	Size	Weight	Colour
922007103	S	35–65 kg (77–143 lb)	Red
922007203	M	55–85 kg (121–187 lb)	Yellow
922007303	L	75–125 kg (165–275 lb)	Green
922007403	XL	120–160 kg (264–352 lb)	Blue
922007504	XXL	150–230 kg (330–506 lb)	Blue



HANDICARE

485 Millway Avenue, Unit 2, Concord, ON, L4K 3V4

Toll Free: 877.304.5438 | Fax: 416.260.0256

customerservice.canada@handicare.com

www.handicare.ca

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

customerservice.northamerica@handicare.com

www.handicareusa.com

North America's largest sling manufacturer. Let us help you find the correct sling for any application.

To learn more, call 1-888-637-8155 or visit www.handicare.ca / handicareusa.com