




Xclusive Straight Stairlift Dimensions



| Perch Seat | Dimensions |
|--|------------|
| A – Minimum swivel radius | 27.5in |
| B – Wall to inner seat back | 3.75in |
| C – Seat depth | 15in |
| D – Minimum folded width to footrest | 13in |
| E – Minimum open width to edge of footrest | 22.125in |
| F – Armrest width (external) | 25.625in |
| G – Armrest width (internal) | 13.75in |
| H – Seat back height | 16.5in |
| I – Footrest to seat height | 16in |
| · One spacer | 18.125in |
| · Two spacers | 19.875in |
| J – Minimum footrest height | 4in |
| Minimum track intrusion into staircase | 7.25in |
| Maximum weight limit | 285lbs |

Some of the above dimensions depend on the characteristics of the staircase in the home. Please consult the sales manual for further information about these dimensions and features of this seat. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.